

# Summary of Papers regarding: Autonomous Tip Conditioning

## Contents

Summary of Papers regarding: Autonomous Tip Conditioning .....	1
Conclusion .....	1
Automated Tip Conditioning for Scanning Tunneling Spectroscopy .....	2
Autonomous In Situ Tip Conditioning via Machine Learning .....	3
Scanbot: An STM Automation Bot .....	4
Automated Scanning Probe Tip State Classification without Machine Learning .....	5
Automated Structure Discovery for Scanning Tunneling Microscopy* .....	6
Bibliography .....	7

## Conclusion

# Automated Tip Conditioning for Scanning Tunneling Spectroscopy

[1] S. Wang, J. Zhu, R. Blackwell, and F. R. Fischer, “Automated Tip Conditioning for Scanning Tunneling Spectroscopy,” *The Journal of Physical Chemistry A*, vol. 125, no. 6, pp. 1384–1390, Feb. 2021, doi: 10.1021/acs.jpca.0c10731.

## Method Overview

Wang et al. (2021) develop a fully automated, closed-loop protocol to prepare atomically sharp scanning-tunneling-microscope (STM) tips at 4.2 K. Their routine alternates between “conditioning” (mechanical pokes) and “assessment” ( $\frac{dI}{dV}$  spectroscopy + machine-learning classification) until two consecutive spectra indicate a high-quality tip, or until all candidate sites are exhausted and the scan area shifts.

## Workflow Steps

### 1. Image Acquisition & Flattening

- Capture a  $100\text{ nm} \times 100\text{ nm}$  topograph at  $V_{\text{bias}} = 50\text{ mV}$ ,  $I_t = 20\text{ pA}$ .
- Remove tilt by fitting a plane to either the entire image (uniform surfaces) or three widely separated flat regions (stepped/molecular surfaces).

### 2. Surface Segmentation

- Build a height histogram.
- Detect peaks and group pixels within  $\pm 0.05\text{ nm}$  of each peak into “terrace” labels.

### 3. Site Selection

- Slide a  $5\text{ nm} \times 5\text{ nm}$  window across the image.
- Select centers of uniformly labeled squares at least  $15\text{ nm}$  apart as conditioning sites.

### 4. Conditioning & Assessment Loop

1. Move the tip to the chosen site; perform a  $2\text{ nm}$  “poke.”
2. Record two  $dI/dV$  spectra (lock-in frequency =  $455\text{ Hz}$ , modulation =  $10\text{ mV}$ ).
3. Normalize the second spectrum over  $-1.5\text{...}2.0\text{ V}$  ( $896$  points).
4. Classify with an AdaBoost model:
  - **Good:** mark success; two consecutive “good”  $\rightarrow$  terminate.
  - **Bad:** repeat poke + spectroscopy at the same site.

### 5. Scan-Area Shift

- If no site yields two good spectra, shift the scan window by  $100\text{ nm}$  in X or Y and repeat from Step 1 until piezo limits are reached.

## Key Findings & Recommendations

- **AdaBoost** offers the best balance of accuracy, precision, recall, and implementation simplicity.
- Deep nets require more computational overhead and careful tuning.
- The closed-loop poke, measure, classify cycle converges to a publication-quality tip in fewer than ten attempts on average.

## Short Method Summary

An automated loop acquires and flattens STM images, segments flat terraces, selects conditioning sites, and iteratively “pokes” and measures  $dI/dV$  spectra at each site. A lightweight AdaBoost classifier judges tip quality; two consecutive “good” calls end conditioning. If all local sites fail, the scan region shifts and the process repeats until a sharp tip is achieved.

# Autonomous In Situ Tip Conditioning via Machine Learning

[2] M. Rashidi and R. A. Wolkow, "Autonomous Scanning Probe Microscopy in Situ Tip Conditioning through Machine Learning," *ACS Nano*, vol. 12, no. 6, pp. 5185–5189, Jun. 2018, doi: 10.1021/acsnano.8b02208.

## Method Overview

Rashidi & Wolkow (2018) present an automated routine for in situ conditioning of STM tips during hydrogen-terminated Si(100) experiments. A convolutional neural network (CNN) analyzes isolated dangling-bond images to detect degraded ("double") tips; upon detection, controlled tip indentations restore sharpness until the CNN confirms a single-atom apex.

## Workflow Steps

### 1. Data Acquisition & Preprocessing

- Collect ~3500 STM sub-images of isolated dangling bonds at  $-1.8\text{ V}$ ,  $50\text{ pA}$  ( $5.6 \times 5.6\text{ nm}^2$ ).
- Resize each to  $28 \times 28\text{ px}$ ; augment by four  $90^\circ$  rotations and mirror ( $\times 8$  total).

### 2. Baseline Benchmark (Pearson Correlation)

- Compute Pearson's coefficient against sharp-tip references.
- Grid-search threshold yields 77 % classification accuracy.

### 3. Model Training & Selection

- Evaluate KNN ( $k = 5$ ), RFC (5000 trees), SVM (RBF kernel,  $C = 500$ ,  $\gamma = 0.5$ ), FCNN (18 layers  $\times$  784 nodes), and CNN.
- Optimal CNN: conv $5 \times 5$  filters ( $30 \rightarrow 40$  channels, stride 1)  $\rightarrow$  max-pool  $2 \times 2 \rightarrow$  dense 128 nodes (ReLU)  $\rightarrow$  softmax output; trained with Adam (lr  $10^{-4}$ ), cross-entropy loss.

### 4. Automated Tip-Conditioning Loop

1. Acquire full-frame STM image; detect and extract dangling-bond patches.
2. Classify each patch via CNN; perform majority voting across  $N$  patches ( $>99\%$  reliability).
3. If tip = "double," perform indentation at a user-preset spot: approach  $700\text{ pm} \rightarrow 1\text{ nm}$  beyond setpoint ( $-1.8\text{ V}$ ,  $50\text{ pA}$ ), incrementing by  $10\text{ pm}$  on failure.
4. Repeat acquisition and classification until CNN outputs "sharp."

### 5. Integration with Atomic Fabrication

- Demonstrated during binary atomic wire patterning: routine paused fabrication only to recondition when tip degraded, then resumed seamlessly.

## Key Findings & Recommendations

- CNN outperforms classical and shallow ML methods, achieving 97 % raw accuracy and  $>99\%$  with voting.
- Majority voting over multiple defects significantly reduces misclassification risk.
- Framework is generalizable to any surface with recurrent atomic features (e.g., defects, adsorbates).

## Short Method Summary

A loop extracts dangling-bond images, classifies tip quality with a CNN, and performs controlled indentations until tip sharpness is confirmed, enabling uninterrupted, autonomous atomic-scale fabrication.

# Scanbot: An STM Automation Bot

[3] J. Ceddia, J. Hellerstedt, B. Lowe, and A. Schiffrin, “Scanbot: An STM Automation Bot,” *Journal of Open Source Software*, vol. 9, no. 99, p. 6028, Jul. 2024, doi: 10.21105/joss.06028.

## Method Overview

Ceddia et al. (2024) introduce Scanbot, a Python-based “robot” that fully automates key STM tasks—tip conditioning, sample surveying, and data acquisition—by coordinating piezoceramic scanners with a real-time camera feed and requiring **Nanonis V5** control software for STM integration.

## Workflow Steps

### 1. DSH Calibration & Tip Tracking

- Initialize camera feed to locate and track the STM tip apex and target positions on both the sample and a clean reference metal.
- Use piezoceramic scanner commands (via Nanonis V5 API) to maneuver the tip between regions.

### 2. Sample Survey & Site Identification

- Acquire a coarse topographic map of the sample area.
- Identify regions of interest (flat, debris-free patches) on both the sample and the clean metal for imaging and conditioning.

### 3. Tip-Shaping Loop

1. **Imprint Generation:** Gently impinge the tip onto the clean reference metal to leave an atomic-scale imprint.
2. **Imprint Imaging:** Scan the imprint region to produce an image that reflects the tip’s geometry.
3. **Quality Assessment:** Measure **area** and **circularity** of the imprint; compare against predefined thresholds.
4. **Conditional Shaping:**
  - If criteria are met → tip deemed “sharp.”
  - If not → perform a more aggressive poke at a new location on the metal and repeat steps 1-4.

### 4. Resumption of Data Acquisition

- Once the imprint satisfies quality metrics for two consecutive assessments, automatically return the tip to the sample of interest and resume STM imaging or spectroscopy.

## Key Findings & Recommendations

1. **Software Compatibility** Scanbot is compatible with any STM system controllable via **Nanonis V5**.
2. **Modular Design** All core functionalities are hook-based, allowing labs to plug in custom imaging or conditioning routines without modifying the Scanbot core.
3. **Imprint-Based Metrics** Quantitative analysis of imprint geometry (area & circularity) provides a robust, microscope-agnostic metric for tip quality.
4. **Performance** Demonstrated reliable convergence to high-quality tips in a handful of shaping cycles, reducing manual overhead.

## Short Method Summary

Scanbot leverages **Nanonis V5** control software to orchestrate a closed-loop sequence, camera-guided tip positioning, imprint-based tip-shaping on a clean reference metal, and quantitative image analysis, automatically restoring and maintaining an atomically sharp STM probe before returning to the sample for uninterrupted, high-resolution data acquisition.

# Automated Scanning Probe Tip State Classification without Machine Learning

[4] D. S. Barker, P. J. Blowey, T. Brown, and A. Sweetman, “Automated Scanning Probe Tip State Classification without Machine Learning,” *ACS Nano*, vol. 18, no. 3, pp. 2384–2394, Jan. 2024, doi: 10.1021/acsnano.3c10597.

## Method Overview

Barker et al. (2024) present a template-matching (TM)–based approach to classify STM tip state directly from a single topographical image without requiring large labeled data sets or machine learning. Their LabVIEW scripts interface with an RC5 **Nanonis** controller to automate image acquisition, classification, and in situ tip conditioning via bias pulses and nano-indentations.

## Workflow Steps

### 1. Image Acquisition & Preprocessing

- Acquire constant-current STM topographs ( $20 \times 20 \text{ nm}^2$ ,  $720 \times 720 \text{ px}$ ) at system-specific biases and setpoints (e.g., 2 V/200 pA for Si(111)  $7 \times 7$ ).
- Flatten images and remove bottom scan lines to avoid creep artifacts.

### 2. Template Matching Classification

#### 1. Cross-Correlation (CC)

- Choose a small reference template (e.g., corner-hole on Si(111)  $7 \times 7$  or dangling bond on B:Si(111)).
- Compute the cross-correlation ratio (CCR) feature map by scanning the template over the input image.
- Aggregate top N CCR peaks to derive a single metric; classify as **Good** or **Bad** based on an empirically determined threshold.

#### 2. Circularity Measurement (for adatoms)

- Locate adatom via CC; normalize and binarize the region at multiple thresholds (0.4-0.7).
- Compute circularity  $C(r) = \frac{\sigma(r)}{\bar{r}}$  across thresholds; classify based on threshold (e.g.,  $C < 0.035$  denotes **Good**).

### 3. Automated Tip Preparation Loop

- Repeatedly acquire images and classify tip state.
- If **Bad**, move 200 nm away, apply bias pulses and indentations of increasing magnitude to condition the tip.
- After a set number of shaping attempts or upon classification as **Good**, return to imaging site or shift scan area using the coarse motor.

## Key Findings & Recommendations

1. **Accuracy & Precision** TM classifier achieves 90% accuracy and >95% true-positive precision (TPP) on Si-based surfaces, comparable to CNN models and human operators.
2. **Minimal Overhead** Requires only a single **Good** reference image; no training or large labeled data sets needed.
3. **Wide Applicability** Effective on systems with repeating surface features; flexible to new substrates by selecting appropriate templates.
4. **Robust Automation** Integrated within LabVIEW and RC5 **Nanonis** for a fully autonomous tip conditioning tool, achieving a **Good** tip in 12 shaping events (10 min).

## Short Method Summary

A LabVIEW–driven tool using an RC5 **Nanonis** controller automates STM tip state classification via cross-correlation and circularity template matching, then performs bias-pulse and indentation conditioning until consecutive **Good** classifications restore an atomically sharp tip for uninterrupted high-resolution imaging.

# Automated Structure Discovery for Scanning Tunneling Microscopy\*

[5] L. Kurki, N. Oinonen, and A. S. Foster, “Automated Structure Discovery for Scanning Tunneling Microscopy,” *ACS Nano*, vol. 18, no. 17, pp. 11130–11138, Apr. 2024, doi: 10.1021/acsnano.3c12654.

## Method Overview

Kurki et al. (2024) present **ASD-STM**, a machine-learning pipeline that predicts atomic structure directly from bond-resolved STM images by training an Attention U-Net on simulated STM/AFM datasets. No real-time STM control is required; the method operates on acquired images

## Workflow Steps

### 1. Synthetic Dataset Generation

- Compute electronic states for ~81,000 organic molecules using FHI-aims (PBE functional).
- Simulate constant-height STM images via PPSTM (Bardeen tunneling theory) and AFM images via PPAFM over a range of tip-sample distances.
- Generate atomic disk descriptors (size  $\propto$  covalent radius, brightness  $\propto$  height) for training

### 2. Machine Learning Model Training

- Formulate as an image-to-image task: STM image  $\rightarrow$  atomic disk descriptor.
- Develop an Attention U-Net with encoder-decoder blocks and attention gating.
- Train in PyTorch for 50 epochs on 180 k training images (with noise and cut-outs), validate on 20 k, and test on 35.5 k images

### 3. Inference & Validation

- Apply the trained model to experimental bond-resolved STM images of various organic molecules.
- Assess performance: ~1 % mean absolute error on simulated images; ~91 % atom-level and ~74 % ring-level accuracy on hydrocarbons.
- Demonstrate qualitative chemical identification (H vs C vs O) and discuss limitations on challenging cases

## Software Requirements

1. FHI-aims for electronic structure calculations
2. PPSTM & PPAFM for STM/AFM simulations
3. PyTorch for ML model implementation and inference

## Key Findings & Recommendations

1. **High Prediction Accuracy** Mean absolute error ~1 % on simulated images; ~91 % atom accuracy on experimental STM
2. **No Real-Time Control Needed** Decouples ML from instrument operation: runs post-acquisition on standard STM images.
3. **Generalizable Pipeline** Synthetic-dataset approach enables extension to diverse molecules without modifying STM control software.
4. **Future Directions** Incorporate simultaneous AFM signals, vary tip-orbital contributions, and tailor datasets for specialized molecular domains.

## Short Method Summary

**ASD-STM** trains an Attention U-Net on a large, simulated STM/AFM image dataset (via FHI-aims, PPSTM, PPAFM), then predicts atomic structure from a single experimental STM image with high accuracy, requiring only PyTorch for inference and no specialized STM control software.

## Bibliography

- [1] S. Wang, J. Zhu, R. Blackwell, and F. R. Fischer, “Automated Tip Conditioning for Scanning Tunneling Spectroscopy,” *The Journal of Physical Chemistry A*, vol. 125, no. 6, pp. 1384–1390, Feb. 2021, doi: 10.1021/acs.jpca.0c10731.
- [2] M. Rashidi and R. A. Wolkow, “Autonomous Scanning Probe Microscopy in Situ Tip Conditioning through Machine Learning,” *ACS Nano*, vol. 12, no. 6, pp. 5185–5189, Jun. 2018, doi: 10.1021/acs.nano.8b02208.
- [3] J. Ceddia, J. Hellerstedt, B. Lowe, and A. Schiffrin, “Scanbot: An STM Automation Bot,” *Journal of Open Source Software*, vol. 9, no. 99, p. 6028, Jul. 2024, doi: 10.21105/joss.06028.
- [4] D. S. Barker, P. J. Blowey, T. Brown, and A. Sweetman, “Automated Scanning Probe Tip State Classification without Machine Learning,” *ACS Nano*, vol. 18, no. 3, pp. 2384–2394, Jan. 2024, doi: 10.1021/acs.nano.3c10597.
- [5] L. Kurki, N. Oinonen, and A. S. Foster, “Automated Structure Discovery for Scanning Tunneling Microscopy,” *ACS Nano*, vol. 18, no. 17, pp. 11130–11138, Apr. 2024, doi: 10.1021/acs.nano.3c12654.